

PLEASE CHECK WWW.MOLEX.COM FOR LATEST PART INFORMATION

Part Number: [19027-0079](#)
Status: **Active**
Description: ATP Perishable Tool Kit used in ATP Presses that Process Product on Tape

Documents:
[RoHS Certificate of Compliance \(PDF\)](#)

General

Product Family	Application Tooling
Series	T1102
Comments	Conductor Anvil (Lower) 19027-0072, Conductor Punch (Upper) 19027-0073, Insulation Anvil (Lower) 19027-0080, Insulation Punch (Upper) 19027-0081
Function	Crimp
Geographic Area	Global
Level of Automation	Semi-Automatic
More Detailed Tech Information	applicationtooling@molex.com
Product Name	N/A
Survey	www.molex.com//review/apptool_survey.html
Tool Type	Tool Kit
UPC	800753041080

Material Info

Old Part Number ATP-BB-475

Reference - Drawing Numbers

Application Tooling Documents
 ATS-190270079, ATS-TOOLING-ID-ATP, TM-640160065, TM-640160065CN, TM-640160065DE, TM-640160065FR, TM-640160065IT, TM-640160065JP, TM-640160065KR, TM-640160065MY, TM-640160065PL, TM-640160065RU, TM-640160065SK, TM-640160065SP, TM-640160065TH, TM-640160065TW, TM-640160065VN




image - Reference only

EU RoHS **China RoHS**

Compliance Status

Not Reviewed

REACH SVHC

Not Reviewed

Low-Halogen Status

Not Reviewed

Need more information on product environmental compliance?

Email productcompliance@molex.com
 For a multiple part number RoHS Certificate of Compliance, [click here](#)

Please visit the [Contact Us](#) section for any non-product compliance questions.

Search Parts in this Series

T1102Series

Use With

ATP-301 press, ATP-101 press, ATP-101W press, ATP-101WDI, ATP-101HD

Application Tooling	FAQ
Description	Product #

VibraKrimp™	0190180013
-------------	------------

Quick Disconnect,
 Female, for 14-16
 (2.00-1.30mm²) Wire,
 Mylar Tape, Tab 4.75
 x 0.51mm

VibraKrimp™	0190180015
-------------	------------

Quick Disconnect,
 Female, for 14-16
 (2.00-1.30mm²) Wire,
 Mylar Tape, Tab 6.35
 x 0.81mm

PLEASE CHECK WWW.MOLEX.COM FOR LATEST PART INFORMATION